AXCENT ACRYLIC™

Repair Denture Base

CADMIUM FREE

Processing Instructions

Self-Curing Denture Resin for Repairs, Relines and Rebases

Complies with ISO 20795-1:2008 Specification Type 2, Class 1
Use with Garreco Axcnt Acrylics Self Cure Acrylic Liquid Monomer

INTENDED USE:

For the fabrication of dental appliances.

CAUTION:

This product is to be used only under the supervision of a licensed dental practitioner.

INDICATION(S) FOR USE:

Garreco Axcnt Acrylics Repair Denture Base is used for the repair of denture bases.

WARNINGS AND PRECAUTIONS:

IF ON SKIN: May be irritating to skin in some sensitive individuals, especially after prolonged and/or repeated contact. IF INHALED: May be irritated by gross overexposure, no matter how generated. IF SWALLOWED: Not expected to cause any harmful effects. Ingestion of large amounts may be irritating if product is swallowed, may cause nausea, headache, vomiting and/or diarrhea. IF IN EYES: May be irritating to the eyes by gross overexposure, no matter how generated. Symptoms of overexposure may include redness, itching, irritation and watering. Keep dust out of eyes.

CONTRAINDICATION(S):

Do not use with patients with a known allergy to methacrylates or other components of the product.

ADVERSE REACTION(S):

Reversible inflammatory changes of the oral mucosa may occur with susceptible persons.

INSTRUCTIONS:

1. Prepare Appliance

Select a shade to match original appliance. Splint the broken appliance using sticky wax to secure pieces. Pour a working cast and allow to set. Coat surface(s) of cast with separating film. With appropriate grinding tools, open the area of fracture to accommodate repair material.

2. Use either the "Salt and Pepper" or the Mixing Jar Technique listed below

a. "Salt and Pepper" Repair Technique

Mount appliance on cast. Apply monomer (as above) to beveled edges, add powder to the moistened parts. Continue to add successive layers of monomer and powder until the repair is "overbuilt". Powder and monomer dispensers are available for this application.

b. Mixing Jar Technique

For a doughy mix, by volume add 2 parts powder to 1 part Self Cure Liquid (6.5 g powder to 5 ml liquid) and spatulate for 20 seconds or until thoroughly mixed. Do not over spatulate. When the mixture has achieved a doughy consistency, it is ready to pack. Wet beveled edges of case with monomer then using an instrument, pack mix to an "overbuilt" state.

3. Curing

Cover the case completely in warm water (115°F/46°C) in a pneumatic pressure pot. Cure for 15-30 minutes at 20 psi depending on thickness of repair. CAUTION: Curing at too high a temperature could result in warpage of the repaired appliance. Bench cool the case to room temperature and finish the overbuilt area to an even surface. Finish and polish as usual.

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Garreco, LLC