

DYNASTY HI-IMPACT POUR TYPE™ DENTURE BASE

Introducing the world's best self cure, high impact denture base.

New from Garreco Dental. Dynasty Hi-Impact Pour Type is the first dental acrylic to achieve high strength and high impact resistance from a self cure formulation, achieving a fracture toughness 35% greater than conventional self cure denture base polymers. It is the only self cure denture base on the market to fulfil the requirements of standard EN ISO 20795-1:2013 for high-impact denture base polymers (Sections 5.2.10 & 5.2.11). This new self cure material allows dental labs to cope with high production while keeping energy costs and labor to a minimum by avoiding the equipment, increased labor, energy use and cure times associated with heat cured materials.



NEW PRODUCT - AVAILABLE THE SECOND QUARTER OF 2017

SOFTLINER RELINE MATERIAL KIT

A self-mixing material that exudes simplicity. The Softliner Reline Base comes in a 2-part cartridge with mixing tips that negates manual measuring and mixing making easy and consistent results. The cartridge base also allows material to never be exposed to air, thus keeping unused product fresh and ready for next application. Long term expiration is appealing to all levels of production.

Indication for use: Relining dentures

Kits Available:

- Softliner Reline Material Kit
- Softliner Application Tools Kit
- Softliner Refill Kit
- Softliner Primer Kit
- Softliner Gun Kit
- Softliner Varnish Kit



PERLADIA POLISHERS DISTRIBUTED IN THE US EXCLUSIVELY BY GARRECO

Perlada Polishers that have been the standard in porcelain polishing for many years and can be used with any porcelain. They are long lasting, generate no heat and eliminates re-glazing. No porcelain chipping and no polishing paste. All PD Polishers are impregnated with diamonds!

Indication for use: Shaping, finishing and polishing with no chipping or heat generation

Items Available:

Knife Edge - 25 mm - Extra Coarse Ochre, Coarse Red, Medium-Coarse Tan, Medium Blue, Fine Green, Extra Fine Purple, High Shine Yellow, Coarse Black

Square Edge - 15 mm - Extra Coarse Ochre, Coarse Red, Medium-Coarse Grey, Medium Blue



SILAFORM™ SILICONE PUTTY

Silaform silicone putty is a one to one mix material which is excellent for a variety of procedures. Silaform 90 (blue) has a Shore A hardness of 90 while Silaform 70 (purple) is more flexible and has a Shore A hardness of 70.

Both materials resist deformation and abrasion normally associated with silicone putty. The absence of expansion and contraction offers endless possibilities for this **economically priced** product from Garreco.

Indication for use: Duplicating material

Sizes Available:

70 Hard - 1kg, 25kg
90 Extra Hard - 1.5kg, 25kg



ACCUSIL™ SILICONE DUPLICATING MATERIAL

A high quality platinum based silicone duplicating material for use with all techniques. Accuracy is constant and predictable at 0.0% expansion/contraction. Excellent strength of 16-18 Shore A hardness. May be poured multiple times with gypsum or investment without deterioration. Suitable for use in dispensers.

Indication for use: Duplicating material

Sizes Available: 1kg, 25kg

