

RESILIENCE™ POLISHING/FINISHING MEDIA

Resilience™ is an all-purpose acrylic polishing media designed for faster and smoother polishing of acrylic appliances. Specially formulated with unique abrasives to cut faster, easier, and smoother than pumice-type products with no scratching or grooving. Resilience saves time and does not require a change in your current technique. Kind to your hands and does not give off a bad odor commonly associated with pumice.

- Easy to use
- Replacement for pumice media
- One grit size does it all
- No odor or bacteria
- No scratching or grooving
- Saves time

Indication for use: Polishing acrylic appliances

Sizes: 2.5lb, 10lb, 20lb



ALUMINUM OXIDE ABRASIVE

Garreco Aluminum Oxide blast abrasives are among the finest in the world and used for divesting and cleaning of appliances. Our 99.6% pure (white) and 96% pure (brown) electro-fused Aluminum Oxide abrasives are machine-sized and computer-scanned to provide uniform particle size. Available in popular sizes and packed with desiccant to prevent moisture absorption.

- Easy to use
- Saves time
- Ultimate performance
- Multiple sizes for different applications

Indication for use: Divesting/cleaning of appliances

Sizes: 2lb, 6 x 2lb, 25lb, 50lb

Particle size (white): 25, 50, 60, 100, 110, 250 microns

Particle size (brown): 50 microns



WALNUT SHELLS FINISHING MEDIA

All-natural Walnut Shells from Garreco are ideal for use after divesting full and partial dentures to easily remove gypsum from the acrylic surface. Available in 20/30 grit only.

- Economical
- Natural

Indication for use: Removing gypsum from acrylic surfaces

Sizes: 20lb, 40lb

Particle size: 20/30 grit



GLASS BEADS ABRASIVE

Glass Beads are used to divest chrome cobalt, non-precious crown and bridge alloys, and all-ceramic cases. May also be used to polish chrome cobalt frameworks. Precision screened and desiccant packaged to prevent moisture.

- Consistent particle size
- Precision screened

Indication for use: Divesting, ceramic cases, polishing chrome cobalt

Sizes: 25lb, 50lb

