



Garreco

ACCUSIL SILICONE DUPLICATING MATERIAL

A high quality platinum based silicone duplicating material for use with all techniques. Accuracy is constant and predictable at 0.0% expansion/contraction. Excellent strength of 16-18 Shore A hardness. May be poured multiple times with gypsum or investment without deterioration. Suitable for use in dispensers.

- Sold in sets or individually
- 0.0% expansion
- High strength

Indication for use: Duplicating material

MIXING RATIO..... 1:1 A/B
 WORKING TIME >5 min
 CURING TIME 30-45 min
 ELONGATION OF TEAR ~310%




ACCUSIL SIZES AVAILABLE

- 1 KG SET (1 ACTIVATOR & 1 BASE) - 6705005
- 25 KG SET (1 ACTIVATOR & 1 BASE) - 6705055

STARTER SET

ACCUSIL STARTER SET - 6700006
 Starter Kit Includes: 1 each 8oz Accusil A & B, 10" Duplicating Ring, Duplicating Cross, Silasavers (1 each: sm, med, lg)



Manufactured by Garreco, LLC
 430 Hiram Road, Heber Springs, AR 72543, USA
 Toll Free: (800) 334-1443
 garreco.com