



Garreco

DIE STONE

APEX RZ 10 TYPE IV-R RESIN REINFORCED DIE STONE

Apex is a resin reinforced natural gypsum die material, available in two expansion variations making it the perfect choice to complement any technique. Perfect for all crown and bridge applications as well as implant or precision attachment cases. Resin additives eliminate brittleness and add superior external hardness. The exclusive Apex flow agent allows for excellent stackability, while promoting free material flow during vibration.

- Exceptionally dense for maximum strength
- High surface abrasion resistance
- Works well with all impression material
- Captures every detail

Indication for use: Crown and bridge; porcelain; implants; cast partials; high-precision cases

SETTING TIME 10-13 min
 WATER/POWDER RATIO 19ml:100g
 COMPRESSIVE STRENGTH 18,500 psi
 EXPANSION 0.10%



	BLUE	GREEN	PINK	IVORY*	WHITE	TAN
120x70g	1301580	1302580	1305080	1303580	1303080	1306580
25LB	1301525	1302525	1305025	1303525	1303025	1306525
50LB	1301550	1302550	1305050	1303550	1303050	1306550

Manufactured by Garreco, LLC
 430 Hiram Road, Heber Springs, AR 72543, USA
 Toll Free: (800) 334-1443
 garreco.com

* **WARNING: Cancer** - www.P65Warnings.ca.gov

Garreco, LLC is an ISO 13485 Facility
 Garreco is a trademark of Garreco, LLC. Heber Springs, AR. ©2019 All rights reserved.