



Garreco

AXCENT ACRYLIC 15 MINUTE DENTURE BASE

15-Minute Denture Base is a revolutionary acrylic material. Simply pack as usual, place into boiling water for 15 minutes, then finish as usual. Cures without porosity and its strength is comparable to most premium products.

- No porosity
- Available in 7 shades
- Heat Cure

Indication for use: Denture base

CURE TYPE 15 Minute Boil
 INJECTABLE No
 MICROWAVABLE No
 POWDER/LIQUID RATIO 3:1



MATCHING LIQUID MONOMER

Liquid Monomer #2 - 15 Minute
 8 OZ BOTTLE - 38000085
 8 OZ CASE (4x8 OZ) - 3800008
 1 QT BOTTLE - 38000305
 1 QT CASE (4x1 QT) - 3800030

Manufactured by Garreco, LLC
 430 Hiram Road, Heber Springs, AR 72543, USA
 Toll Free: (800) 334-1443
 garreco.com

HEAT CURE



DENTURE BASE COLORS AND SIZES AVAILABLE

PINK	LIGHT FIBER	DARK FIBER
1 LB JAR - 37080015 1 LB CASE (4x 1LB) - 3708001 5 LB JAR - 37088055 5 LB CASE (4x 5LB) - 3708805 20 LB PAIL - 3708820	1 LB JAR - 37302015 1 LB CASE (4x1 LB) - 3730201 5 LB JAR - 37302055 5 LB CASE (4x5 LB) - 3730205 20 LB PAIL - 3730220	1 LB JAR - 37003015 1 LB CASE (4x1 LB) - 3700301 5 LB JAR - 37003055 5 LB CASE (4x5 LB) - 3700305 20LB PAIL - 3700320
LIGHT CHARACTERIZED	PURPLE FIBER	SHADE 99
1 LB JAR - 37304015 1 LB CASE (4x 1 LB) - 3730401 5 LB JAR - 37304055 5 LB CASE (4x5 LB) - 3730405 20LB PAIL - 3730420	1 LB JAR - 37334015 1 LB CASE (4x1 LB) - 3733401 5 LB JAR - 37334055 5 LB CASE (4x5 LB) - 3733405 20LB PAIL - 3733420	1 LB JAR - 37388015 1 LB CASE (4x1 LB) - 3738801 5 LB JAR - 37388055 5 LB CASE (4x5 LB) - 3738805 20LB PAIL - 3738820
CLEAR		
1 LB JAR - 37301015 1 LB CASE (4x1 LB) - 3730101 5 LB JAR - 37301055 5 LB CASE (4x5 LB) - 3730105 20 LB PAIL - 3730120		