



Garreco

AXCENT ACRYLIC REPAIR DENTURE BASE

Repair Acrylic is color matched to all Garreco denture base materials. Contains special additives to prevent "yellowing" and may be used for any acrylic repair.

- Works with all acrylics
- Available in 7 shades
- Non-yellowing

Indication for use: Denture base repair

CURE TYPE Self
 INJECTABLE Yes
 MICROWAVABLE Yes
 POWDER/LIQUID RATIO 3:1



MATCHING LIQUID MONOMER

Liquid Monomer #1 - Heat Cure
 8 OZ BOTTLE - 3780008S
 8 OZ CASE (4x8 OZ) - 3780008
 1 QT BOTTLE - 3780030S
 1 QT CASE (4x1 QT) - 3780030

Manufactured by Garreco, LLC
 430 Hiram Road, Heber Springs, AR 72543, USA
 Toll Free: (800) 334-1443
 garreco.com

HEAT CURE



DENTURE BASE COLORS AND SIZES AVAILABLE

PINK	LIGHT FIBER	DARK FIBER
1 LB JAR - 3018001S 1 LB CASE (4x 1LB) - 3018001 5 LB JAR - 3018005S 5 LB CASE (4x 5LB) - 3018005 20 LB PAIL - 3018020	1 LB JAR - 3010201S 1 LB CASE (4x1 LB) - 3010201 5 LB JAR - 3010205S 5 LB CASE (4x5 LB) - 3010205 20 LB PAIL - 3010220	1 LB JAR - 3010301S 1 LB CASE (4x1 LB) - 3010301 5 LB JAR - 3010305S 5 LB CASE (4x5 LB) - 3010305 20LB PAIL - 3010320
LIGHT CHARACTERIZED	PURPLE FIBER	SHADE 99
1 LB JAR - 3010401S 1 LB CASE (4x 1 LB) - 3010401 5 LB JAR - 3010405S 5 LB CASE (4x5 LB) - 3010405 20LB PAIL - 3010420	1 LB JAR - 3013401S 1 LB CASE (4x1 LB) - 3013401 5 LB JAR - 3013405S 5 LB CASE (4x5 LB) - 3013405 20LB PAIL - 3013420	1 LB JAR - 3018801S 1 LB CASE (4x1 LB) - 3018801 5 LB JAR - 3018805S 5 LB CASE (4x5 LB) - 3018805 20LB PAIL - 3018820
CLEAR		
1 LB JAR - 3010101S 1 LB CASE (4x1 LB) - 3010101 5 LB JAR - 3010105S 5 LB CASE (4x5 LB) - 3010105 20 LB PAIL - 3010120		