



Garreco

AXCENT ACRYLIC STANDARD DENTURE BASE

A new manufacturing process enables Garreco to offer this impact resistant material at substantial cost savings over lower quality competitors.

- Impact resistant
- Available in 7 shades
- Heat Cure

Indication for use: Denture base

CURE TYPE	Heat
INJECTABLE	Yes
MICROWAVABLE	Yes
POWDER/LIQUID RATIO	3:1



MATCHING LIQUID MONOMER

- Liquid Monomer #1 - Heat Cure**
 8 OZ BOTTLE - 37800085
 8 OZ CASE (4x8 OZ) - 3780008
 1 QT BOTTLE - 37800305
 1 QT CASE (4x1 QT) - 3780030

Manufactured by Garreco, LLC
430 Hiram Road, Heber Springs, AR 72543, USA
Toll Free: (800) 334-1443
garreco.com

HEAT CURE



DENTURE BASE COLORS AND SIZES AVAILABLE

PINK	LIGHT FIBER	DARK FIBER
1 LB JAR - 37080015	1 LB JAR - 37002015	1 LB JAR - 37003015
1 LB CASE (4x 1LB) - 3708001	1 LB CASE (4x1 LB) - 3700201	1 LB CASE (4x1 LB) - 3700301
5 LB JAR - 37088055	5 LB JAR - 37002055	5 LB JAR - 37003055
5 LB CASE (4x5 LB) - 3708805	5 LB CASE (4x5 LB) - 3700205	5 LB CASE (4x5 LB) - 3700305
20 LB PAIL - 3708820	20 LB PAIL - 3700220	20 LB PAIL - 3700320
LIGHT CHARACTERIZED	PURPLE FIBER	SHADE 99
1 LB JAR - 37004015	1 LB JAR - 37034015	1 LB JAR - 37088015
1 LB CASE (4x 1 LB) - 3700401	1 LB CASE (4x1 LB) - 3703401	1 LB CASE (4x1 LB) - 3708801
5 LB JAR - 37004055	5 LB JAR - 37034055	5 LB JAR - 37088055
5 LB CASE (4x5 LB) - 3700405	5 LB CASE (4x5 LB) - 3703405	5 LB CASE (4x5 LB) - 3708805
20LB PAIL - 3700420	20LB PAIL - 3703420	20LB PAIL - 3708820
CLEAR		
1 LB JAR - 37001015		
1 LB CASE (4x1 LB) - 3700101		
5 LB JAR - 37001055		
5 LB CASE (4x5 LB) - 3700105		
20 LB PAIL - 3700120		