



Garreco

HEAT CURE

DYNASTY ACRYLIC HI-IMPACT DENTURE BASE

Dynasty Hi-Impact is a cadmium free denture base acrylic with ultimate strength and flexibility. It is dimensionally accurate, color stable, requires no change in technique, and has excellent working characteristics. Dynasty Hi-Impact polishes to a dense, high gloss finish.

- Cadmium free
- High impact resistance

Indication for use: Denture base

CURE TYPE Heat
 INJECTABLE Yes
 MICROWAVABLE Yes
 POWDER/LIQUID RATIO 3:1



DENTURE BASE COLOR AND SIZES AVAILABLE

PINK

- 1 LB JAR - 37488015
- 1 LB CASE (4x 1LB) - 3748801
- 5 LB JAR - 37488055
- 5 LB CASE (4x 5LB) - 3748805
- 20 LB PAIL - 3748820



MATCHING LIQUID MONOMER

- Liquid Monomer #1 - Heat Cure**
- 8 OZ BOTTLE - 37800085
 - 8 OZ CASE (4x8 OZ) - 3780008
 - 1 QT BOTTLE - 37800305
 - 1 QT CASE (4x1 QT) - 3780030

Manufactured by Garreco, LLC
 430 Hiram Road, Heber Springs, AR 72543, USA
 Toll Free: (800) 334-1443
 garreco.com