



# Garreco

## SELF CURE

### DYNASTY ACRYLIC HI-IMPACT POUR TYPE DENTURE BASE

Dynasty Hi-Impact Pour Type is the first dental acrylic to achieve high strength and high impact resistance from a self cure formulation, achieving a fracture toughness 35% greater than conventional self cure denture base polymers. This new self cure material allows dental labs to cope with high production while keeping energy costs and labor to a minimum by avoiding the equipment, increased labor, energy use and cure times associated with heat cured materials.

- High strength
- High impact resistance

**Indication for use:** Denture base

CURE TYPE ..... Self  
 INJECTABLE ..... Yes  
 MICROWAVABLE ..... No  
 POWDER/LIQUID RATIO ..... 2:1



#### DENTURE BASE COLOR AND SIZES AVAILABLE

- PINK**
- 
- 1 LB JAR - 39388015
  - 1 LB CASE (4x 1LB) - 3938801
  - 5 LB JAR - 39388055
  - 5 LB CASE (4x 5LB) - 3938805
  - 20 LB PAIL - 3938820



#### MATCHING LIQUID MONOMER Liquid Monomer for Hi-Impact Pour Type Acrylic

- 4 OZ BOTTLE - 3940004S
- 4 OZ CASE (4x4 OZ) - 3940004
- 8 OZ BOTTLE - 3940008S
- 8 OZ CASE (4x8 OZ) - 3940008
- 1 QT BOTTLE - 3940030S
- 1 QT CASE (4x1 QT) - 3940030

Manufactured by Garreco, LLC  
430 Hiram Road, Heber Springs, AR 72543, USA  
Toll Free: (800) 334-1443  
garreco.com