



Garreco

BASE STONE

EXPRESS SET TYPE IV SPECIALTY STONE

ExpressSet is a Type IV specialty stone with an extremely fast set time, smooth mixing consistency, and uniform color matrix. Models can be released from the impression in 5 minutes, producing a precision cast expressly suited for mouth guards, custom trays, study casts, models for opto-electronic/ laser scanning with powder, temporaries and denture repairs.

- High compressive strength
- Ample working time
- Moderate expansion (0.15%)

Indication for use: Crown and bridge; porcelain; implants; cast partials; high-precision cases

SETTING TIME 60-90 sec.
 WATER/POWDER RATIO 23ml:100g
 COMPRESSIVE STRENGTH 14,000 psi
 EXPANSION 0.15%



BUFF



PINK



80x140g	1281092	1285092
25LB	1281025	1285025
50LB	1281050	1285050

Manufactured by Garreco, LLC
 430 Hiram Road, Heber Springs, AR 72543, USA
 Toll Free: (800) 334-1443
 garreco.com