



# Garreco

## BASE STONE

### FASTONE TYPE III DENTAL BASE STONE

This dental stone offers the same unique characteristics as our Tecstone product, with a faster setting time. Used for full or partial denture models, denture flasking procedures, and recommended for use in basing crown and bridge models.

- All-purpose model stone
- Smooth, strong casts
- Low expansion (0.12%) for accuracy

**Indication for use:** Full or partial denture models; flasking; basing crown and bridge models

SETTING TIME ..... 10-12 min  
 WATER/POWDER RATIO ..... 30ml:100g  
 COMPRESSIVE STRENGTH ..... 9,000 psi  
 EXPANSION ..... 0.12%



BLUE



25LB 1101525  
50LB 1101550

BUFF



1101025  
1101050

WHITE



1103025  
1103050

GOLDEN



1106725  
1106750

Manufactured by Garreco, LLC  
430 Hiram Road, Heber Springs, AR 72543, USA  
Toll Free: (800) 334-1443  
garreco.com