



Garreco

GLASS BEADS ABRASIVE MATERIAL

Glass Beads are used to divest chrome cobalt, non-precious crown and bridge alloys, and all-ceramic cases. May also be used to polish chrome cobalt frameworks. Precision screened and desiccant packaged to prevent moisture.

—Consistent particle size

—Precision screened

Indication for use: Divesting, ceramic cases, polishing chrome cobalt

Particle size: 50, 80 microns



GLASS BEADS TYPES AND SIZES AVAILABLE

50 MICRON

25 LB CARTON - 2302025

25 LB PAIL - 2352025

50 LB CARTON - 2302050

50 LB PAIL - 2352050

80 MICRON

25 LB CARTON - 2303525

25 LB PAIL - 2353525

50 LB CARTON - 2303550

50 LB PAIL - 2353550

Manufactured by Garreco, LLC
430 Hiram Road, Heber Springs, AR 72543, USA
Toll Free: (800) 334-1443
garreco.com