



Garreco

BASE STONE

LABSTONE TYPE III DENTAL BASE STONE

The "workhorse" of the gypsum family. Designed for ultimate versatility, Labstone is recommended for use in models and flasking, in the production of full or partial dentures, and in basing crown and bridge models. Creamy mix, moderate compressive strength, and superior outward hardness are all characteristics of this unique formula.

- All-purpose model stone for extremely smooth, strong casts
- 8,500 psi compressive strength offers quality models
- Easy trimming saves time
- Unique particle size offers consistent mixing and setting time

Indication for use: Full or partial dentures; all-purpose model stone; flasking; counter models; basing crown and bridge models

SETTING TIME 11-14 min
 WATER/POWDER RATIO 30ml:100g
 COMPRESSIVE STRENGTH 8,500 psi
 EXPANSION 0.12%



BLUE



25LB 1051525
50LB 1051525

BUFF



1051025
1051050

WHITE



1053025
1053050

Manufactured by Garreco, LLC
430 Hiram Road, Heber Springs, AR 72543, USA
Toll Free: (800) 334-1443
garreco.com