



Garreco

DIE STONE

NATURAL ROCK EH TYPE IV-R RESIN REINFORCED DIE STONE

NaturalRock EH is a traditional, creamy, smooth, resin reinforced die stone which exhibits high compressive strength and very low expansion.

- Exceptionally dense for maximum strength
- Very low expansion (0.08%)

Indication for use: Crown and bridge; porcelain; implants; cast partials; high-precision cases

SETTING TIME	10-13 min
WATER/POWDER RATIO	20ml:100g
COMPRESSIVE STRENGTH	18,000 psi
EXPANSION	0.08%



LT. TAN



120x70g	1446780
25LB	1446725
50LB	1446750

Manufactured by Garreco, LLC
430 Hiram Road, Heber Springs, AR 72543, USA
Toll Free: (800) 334-1443
garreco.com