



Garreco

ORTHODONTIC

ORTHOSTONE TYPE III BASE STONE

OrthoStone exhibits the strength and accuracy of a quality base stone, plus special additives resulting in a brilliant white color. The superior color, working properties, and hardness make this the product of choice for those specializing in orthodontics.

- Super White color
- Ideal hardness

Indication for use: Orthodontic procedures

SETTING TIME	13-15 min.
WATER/POWDER RATIO	30ml:100g
COMPRESSIVE STRENGTH	9,000 psi
EXPANSION	0.12%



WHITE*



25LB	1253025
50LB	1253050

Manufactured by Garreco, LLC
 430 Hiram Road, Heber Springs, AR 72543, USA
 Toll Free: (800) 334-1443
 garreco.com

* **WARNING: Cancer** - www.P65Warnings.ca.gov

Garreco, LLC is an ISO 13485 Facility
 Garreco is a trademark of Garreco, LLC. Heber Springs, AR. ©2019 All rights reserved.