



# Garreco

## OPERATORY GYPSUM

### PRISTINE® III TYPE III BASE STONE

With a broad range of applications this Type III stone may be used for full and partial denture models, denture flasking procedures, and for basing crown and bridge models. Has a fresh, clean scent to mask odors picked up by impression materials which can be transferred to working models. Comes in stackable pails with flip top re-closable lid and a convenient scoop.

- Fast set time
- Fresh, clean scent
- Scoop in each container
- Moderate expansion (0.12%)

**Indication for use:** Gypsum designed for dental office applications

SETTING TIME ..... 8-10 min  
 WATER/POWDER RATIO ..... 30ml:100g  
 COMPRESSIVE STRENGTH ..... 9,000 psi  
 EXPANSION ..... 0.12%



GOLDEN



1026725

WHITE



1023025

25LB PAIL

Manufactured by Garreco, LLC  
430 Hiram Road, Heber Springs, AR 72543, USA  
Toll Free: (800) 334-1443  
garreco.com