



Garreco

OPERATORY GYPSUM


PRISTINE® V TYPE V DIE STONE

Versatility for all your fixed and removable models. This Type V die stone resists normal chipping and abrasion. Has a fresh, clean scent to mask odors picked up by impression materials which can be transferred to working models. Comes in stackable pails with flip top re-closable lid and a convenient scoop.

- Fast set time
- Fresh, clean scent
- Scoop in each container
- Moderate expansion (0.12%)

Indication for use: Gypsum designed for dental office applications

SETTING TIME 10-13 min
 WATER/POWDER RATIO 21ml:100g
 COMPRESSIVE STRENGTH 15,000 psi
 EXPANSION 0.15%



25LB PAIL 1382525 1385025

Manufactured by Garreco, LLC
 430 Hiram Road, Heber Springs, AR 72543, USA
 Toll Free: (800) 334-1443
 garreco.com