



# Garreco

## DIE STONE

### SATINSTONE TYPE IV DIE STONE

Comparable to Whip Mix's Silky-Rock, SatinStone™ is distinguished by exceptional toughness and an extraordinary "satin-smooth" finish that works well with all types of impression materials. Introducing only the finest raw materials allows this product to have a creamy and effortless mix each and every time. SatinStone™ also has a low expansion which makes it one of the most the most unique die stones on the market.

- Exceptionally tough maximum strength
- Creamy, satin-smooth finish
- Low expansion (0.09%-0.11%) means greater accuracy

**Indication for use:** Crown and bridge; porcelain; implants; cast partials; high-precision cases

SETTING TIME ..... ~10 min.  
 WATER/POWDER RATIO ..... 23ml:100g  
 COMPRESSIVE STRENGTH ..... 15,000 psi  
 EXPANSION ..... 0.09%-0.11%



ROSE



WHITE



25LB	1175025	1173025
50LB	1175050	1175050

Manufactured by Garreco, LLC  
 430 Hiram Road, Heber Springs, AR 72543, USA  
 Toll Free: (800) 334-1443  
 garreco.com