



Garreco

BASE STONE

TECSTONE TYPE III DENTAL BASE STONE

Widely recognized as one of the best stones on the market and designed for ultimate versatility, Tecstone is recommended for use in models and flasking, in the production of full or partial dentures, and in basing crown and bridge models. A creamy mix, moderate compressive strength, and superior outward hardness are all characteristics of this unique formula.

- All-purpose model stone
- Smooth, strong casts
- Low expansion (0.12%) for accuracy

Indication for use: Full or partial denture models; flasking; basing crown and bridge models

SETTING TIME	13-15 min
WATER/POWDER RATIO	30ml:100g
COMPRESSIVE STRENGTH	9,000 psi
EXPANSION	0.12%



BLUE



25LB 1011525
50LB 1011550

BUFF



1011025
1011050

WHITE



1013025
1013050

Manufactured by Garreco, LLC
430 Hiram Road, Heber Springs, AR 72543, USA
Toll Free: (800) 334-1443
garreco.com