



# Garreco

## ULTRAVEST UNIVERSAL INVESTMENT MATERIAL

UltraVest excels when utilized with any crown and bridge alloy (except titanium) or any ceramic system. This phosphate bonded, graphite-free formula has a micro-fine particle matrix which allows precise attention to detail. This unique product may be used with either conventional or fast casting techniques and is designed for use with ring or ringless systems. UltraVest may be the only investment you will ever need!

- Carbon-free
- Ring or ringless castings
- All alloys (except titanium)
- Conventional or fast casting
- Liquid/powder sold separately

**Indication for use:** Investment for all crown and bridge alloys

Manufactured by Garreco, LLC  
430 Hiram Road, Heber Springs, AR 72543, USA  
Toll Free: (800) 334-1443  
garreco.com



### ULTRAVEST SIZES AVAILABLE

- 50x100 G POUCHES - 6100012
- 100x90 G POUCHES - 6109023




### MATCHING INVESTMENT LIQUID

- 1 L BOTTLE - 6200003S
- 1 L CASE (4x1 L) - 6200412
- 1 L CASE (16x1 L) - 6201647

Requires 1 liter of expansion liquid per 5 kilograms of investment material

### LIQUID REQUIREMENTS

- 50x100 G - 1 LITER
- 100x90 G - 2 LITERS

\*  **WARNING: Cancer** - [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov)

Garreco, LLC is an ISO 13485 Facility  
Garreco is a trademark of Garreco, LLC. Heber Springs, AR. ©2019 All rights reserved.