



Garreco

VITEX INVESTMENT MATERIAL

Vitex (Ethyl Silicate) partial denture framework investment material has been developed especially for producing chrome cobalt castings. This product offers precise details due to an ultra fine particle size and unique two part blend of liquids. May be used for all high heat techniques. Formulated to mix and pour easily, Vitex investment maintains excellent integrity when used in conjunction with Vitex Model Hardener.

- Precise detail
- Easy mix/pour
- All high heat techniques

Indication for use: Investment material



VITEX INVESTMENT SIZES AVAILABLE

- | | |
|---|--|
| 100 LB KIT | 50 LB KIT |
| ⊘ 100 LB KIT (3 CARTONS) - 5007000
Includes: 2 x 50lb Investment Sand,
2 x 118oz "A" Liquid, and 2 x 14.8oz
"B" Liquid | 50 LB KIT (2 CARTONS) - 5006000
Includes: 1 x 50lb Investment Sand,
4 x 28oz "A" Liquid, and 4 x 3.7oz "B"
Liquid |
| ⊘ 100 LB KIT LIQUID ONLY - 5007000
Includes: 2 x 118oz "A" Liquid, and 2
x 14.8oz "B" Liquid | VITEX LIQUID |
| VITEX SAND | ⊘ 118OZ "A" LIQUID - 5009000 |
| 50LB SAND - 5008000 | 14.8OZ "B" LIQUID - 5011000 |
| | 28OZ "A" LIQUID ONLY - 5010000 |
| | 3.7OZ "B" LIQUID ONLY - 5012000 |

⊘ Denotes an item that cannot be shipped by FedEx Ground; must be shipped FedEx Freight only



HAZMAT ITEM
Special shipping conditions
may apply to these items.

Manufactured by Garreco, LLC
430 Hiram Road, Heber Springs, AR 72543, USA
Toll Free: (800) 334-1443
garreco.com

* **WARNING: Cancer** - www.P65Warnings.ca.gov

Garreco, LLC is an ISO 13485 Facility
Garreco is a trademark of Garreco, LLC. Heber Springs, AR. ©2019 All rights reserved.