



PLASTER

LAB PLASTER
TYPE II PLASTER

A fine white plaster with low setting expansion for models, casts, flasking and standard lab use. Available in Lab Plaster Regular or Lab Plaster Fast Set.

- Low expansion
- Fine, smooth material

Indication for use: Models, casts, flasking and standard lab use

Standard Sizes Available: 25lbs, 50lbs
Additional sizes available for custom orders of 2,000 pounds or more. Contact us for more details.

SETTING TIME	6-8 min (Regular)
	4-6 min (Fast Set)
WATER/POWDER RATIO	47ml:100g
COMPRESSIVE STRENGTH	3,500 psi
EXPANSION	0.18%



REGULAR FAST



25LB	1500025PM	1580025PM
50LB	1500050PM	1580050PM

Manufactured by Garreco, LLC
430 Hiram Road, Heber Springs, AR 72543, USA
Toll Free: (800) 334-1443
garreco.com