



ORTHODONTIC

ORTHO PLASTER
TYPE III BASE STONE

A fine pure white stone recommended for orthodontic casts and working models.

- Perfect for orthodontic cases
- Pure white color

Indication for use: Orthodontic casts and working models

Standard Sizes Available: 25lbs, 50lbs
Additional sizes available for custom orders of 2,000 pounds or more. Contact us for more details.

SETTING TIME	12-14 min
WATER/POWDER RATIO	37ml:100g
COMPRESSIVE STRENGTH	4,800 psi
EXPANSION	0.19%



WHITE



- 25LB 1600025PM
- 50LB 1600050PM

Manufactured by Garreco, LLC
430 Hiram Road, Heber Springs, AR 72543, USA
Toll Free: (800) 334-1443
garreco.com

WARNING: Cancer - www.P65Warnings.ca.gov

Garreco, LLC is an ISO 13485 Facility
Garreco and Pemaco are trademarks of Garreco, LLC. Heber Springs, AR. ©2019 All rights reserved.