



Garreco

PLASTER

ARTICULATING PLASTER TYPE III PLASTER

Articulating Plaster exhibits similar characteristics to our Dental Plasters, but with a greatly reduced set time. Recommended for articulation use.

- Smooth mixing
- Low expansion (0.19%)
- Excellent working characteristics

Indication for use: Articulation

SETTING TIME	3-5 min.
WATER/POWDER RATIO	45ml:100g
COMPRESSIVE STRENGTH	3,500 psi
EXPANSION	0.19%



WHITE



25LB	1560025
50LB	1560050

Manufactured by Garreco, LLC
 430 Hiram Road, Heber Springs, AR 72543, USA
 Toll Free: (800) 334-1443
 garreco.com