



CROWN & BRIDGE



VEST-IT

Vest-it is an asbestos-free, high heat refractory material used to create custom firing pegs.

Indication for use: Creating custom firing pegs

1 SYRINGE (10CC) - 6840002




HONEYCOMB TRAYS

Specifically designed high heat firing trays for use with all furnaces and techniques. Available with ceramic or metal pins..

Indication for use: Platform for firing crowns and bridges

2 TRAYS WITH CERAMIC PINS - 7760801

2 TRAYS WITH METAL PINS - 7760802



FLEX TRAY

This non-stick firing tray combines the best features of rigid trays and soft liners. You'll love the flexibility, durability, and multiple uses for this tray

Indication for use: Flexible tray material for firing crowns and bridges

1 BAG (4 TRAYS/BAG) - 6840001



DISPOSABLE PLUNGERS

Universal disposable plungers for pressing ceramics. Pre-heat not required. Available in 2 sizes (2 gram and 3 gram).

Indication for use: Pressing ceramics

2 G CARTON (50 PCS) - 6080000

3 G CARTON (50 PCS) - 6080300



CRUCIBLES

Type A, non-contaminating, quartz crucibles. Available in slotted or hooded. Sold in pairs.

Indication for use: Melting alloys

HOODED, 2/PACK - CRU-H

SLOTTED, 2/PACK - CRU-S