



Garreco

DEBUBBLIZER SURFACTANT

A superior alcohol-free surfactant for reducing surface tension on wax patterns, impressions, and models.

- Alcohol-free
- Reduces surface tension

Indication for use: Reducing surface tension on wax patterns, impressions, and models



DEBUBBLIZER SIZES AVAILABLE

- 8 OZ BOTTLE - 35515085
- 8 OZ CASE (4x8 OZ) - 3551508
- 1 QT BOTTLE - 35515305
- 1 QT CASE (4x1 QT) - 3551530

Manufactured by Garreco, LLC
430 Hiram Road, Heber Springs, AR 72543, USA
Toll Free: (800) 334-1443
garreco.com