



EXCALIBUR

TYPE IV HIGH-STRENGTH DIE STONE

With its smooth mixing characteristics and ample setting time, Excalibur produces a “textbook” model every time. Die trimming becomes a worry-free procedure due to the ultra-fine gypsum base, which prevents nicks and abrasions caused by brittleness. Recommended for use as crown and bridge or partial denture models. Suitable for use with all impression materials. Gypsum hardeners are not recommended.

- Exceptionally dense for maximum strength
- High surface abrasion resistance
- Works well with all impression material
- Captures every detail

Indication for use: Crown and bridge; porcelain; implants; cast partials; high-precision cases

SETTING TIME 12-15 min
 WATER/POWDER RATIO 22ml:100g
 COMPRESSIVE STRENGTH 18,000 psi
 EXPANSION 0.09%

	BLUE	GREEN	VIOLET	GREY	WHITE	TAN
120x70g	1351580	1352580	1352080	1354580	1353080	1356580
25LB	1351525	1352525	1352025	1354525	1353025	1356525
50LB	1351550	1352550	1352050	1354550	1353050	1356550

Manufactured by Garreco, LLC
 430 Hiram Road, Heber Springs, AR 72543, USA
 Toll Free: (800) 334-1443
garreco.com



DIE STONE