



## BASE STONE

### FLASK IT

#### TYPE III DENTAL BASE STONE

The newest member of the base stone family, allows the technician to make the 2nd and 3rd pours during flasking without vibration. The result is an exceptionally smooth surface which greatly reduces air pockets, bubbles, and valuable time spent removing gypsum around teeth and under-cut areas. This is a must-try for all denture labs.

- Smooth mixing
- Low expansion (0.09%)
- Excellent working characteristics

**Indication for use:** Flasking procedures

SETTING TIME ..... 9-12 min  
WATER/POWDER RATIO ..... 33ml:100g  
COMPRESSIVE STRENGTH ..... 9,000 psi  
EXPANSION ..... 0.09%

PINK



25LB 1235025  
50LB 1235050



Manufactured by Garreco, LLC  
430 Hiram Road, Heber Springs, AR 72543, USA  
Toll Free: (800) 334-1443  
garreco.com