



Garreco

BASE STONE

FLASKING STONE TYPE III DENTAL BASE STONE

A moderate strength, low expansion stone, formulated for denture flasking procedures. Allows rapid heat transfer during boil-out, while maintaining optimum strength during packing. Divests quickly and easily.

- Divests quickly
- Rapid heat transfer

Indication for use: Basing crown and bridge models

SETTING TIME 15-18 min
 WATER/POWDER RATIO 35ml:100g
 COMPRESSIVE STRENGTH 6,500 psi
 EXPANSION 0.12%



WHITE



25LB 1203025
50LB 1203050

BUFF



1201025
1201050

Manufactured by Garreco, LLC
430 Hiram Road, Heber Springs, AR 72543, USA
Toll Free: (800) 334-1443
garreco.com