



Garreco

BASE STONE

MOUNTING STONE TYPE III DENTAL BASE STONE

A high strength, fast setting, low expansion, articulating stone. Originally developed for implant cases where a low expansion articulating material is desired, Many professionals choose this product as a rapid setting alternative material for denture models and crown and bridge bases.

- All-purpose model stone
- Smooth, strong casts
- Low expansion (0.09%) for accuracy

Indication for use: Full or partial denture models; flasking; basing crown and bridge models

SETTING TIME 3-5 min
 WATER/POWDER RATIO 30ml:100g
 COMPRESSIVE STRENGTH 9,000 psi
 EXPANSION 0.09%



WHITE



25LB 1153025
50LB 1153050

GREY



1154525
1154550

Manufactured by Garreco, LLC
430 Hiram Road, Heber Springs, AR 72543, USA
Toll Free: (800) 334-1443
garreco.com