



## NATURAL AGATE

### PALETTE AND SPATULA FOR PORCELAIN

Polished agates provide an ideal platform for mixing dental porcelain. Their varicolor bands add beauty and uniqueness to these one of a kind natural palettes. Agates provide a calming influence, improving perception and concentration, thus enhancing the ceramist's creativity. Agates offer years of use because they do not abrade like porcelain or glass trays/slabs.

The Agate Spatula is an agate instrument ideal for mixing porcelain. The spatulas are metal-free to prevent metal from contaminating the porcelain, plus both ends can be used for mixing making it the perfect accompaniment to a natural agate palette.

- Pleasing to the eye
- Abrasion proof
- Non-contaminating
- One of a kind
- Years of service


**Indication for use:** Tray and spatula for mixing dental porcelain

Manufactured by Garreco, LLC  
430 Hiram Road, Heber Springs, AR 72543, USA  
Toll Free: (800) 334-1443  
garreco.com



#### AGATE ITEMS AVAILABLE

AGATE PALETTE/SPATULA COMBO - 6880010  
AGATE MIXING PALETTE - 680008  
AGATE SPATULA - 6850009



#### HUNDREDS OF AGATES TO CHOOSE FROM!

To view our Natural Agate Gallery, visit us online at [Garreco.com/products/crown-bridge/agate/](http://Garreco.com/products/crown-bridge/agate/) or scan this QR code now