



Garreco

ORTHODONTIC

ORTHOPLASTER TYPE III ORTHODONTIC PLASTER

High strength plaster designed for orthodontic applications. Orthoplaster combines smooth mixing, low expansion and ample setting time. Contains special whitening additives for a brilliant white color.

- Super White color
- Regular or fast set

Indication for use: Orthodontic appliances

SETTING TIME 10-13 min. (Reg); 6-8 min. (fast)
 WATER/POWDER RATIO 35ml:100g
 COMPRESSIVE STRENGTH 6,500 psi
 EXPANSION 0.12%




WHITE*



25LB REG	1600025
50LB REG	1600050
25LB FAST	1620025
50LB FAST	1620050

Manufactured by Garreco, LLC
 430 Hiram Road, Heber Springs, AR 72543, USA
 Toll Free: (800) 334-1443
 garreco.com

*  **WARNING: Cancer** - www.P65Warnings.ca.gov

Garreco, LLC is an ISO 13485 Facility
 Garreco is a trademark of Garreco, LLC. Heber Springs, AR. ©2019 All rights reserved.