



Garreco

OPERATORY GYPSUM

PRISTINE® II TYPE II PLASTER

A fast setting Type II model and articulating plaster exhibiting excellent working characteristics. Has a fresh, clean scent to mask odors picked up by impression materials which can be transferred to working models. Comes in stackable pails with flip top re-closable lid and a convenient scoop.

- Fast set time
- Fresh, clean scent
- Scoop in each container
- Moderate expansion (0.19%)

Indication for use: Gypsum designed for dental office applications

| | |
|----------------------------|-----------|
| SETTING TIME | 6-10 min |
| WATER/POWDER RATIO | 45ml:100g |
| COMPRESSIVE STRENGTH | 3,500 psi |
| EXPANSION | 0.19% |



WHITE



20LB PAIL 1500020

Manufactured by Garreco, LLC
430 Hiram Road, Heber Springs, AR 72543, USA
Toll Free: (800) 334-1443
garreco.com