



Garreco

BASE STONE


TECSTONE FL - FLOWABLE TYPE III DENTAL BASE STONE

This unique new product from Garreco is specifically designed for basing crown & bridge models. Tecstone FL – Flowable pours into base formers without vibration and yields a bubble-free base. Tecstone FL – Flowable far exceeds the strength of conventional basestones while offering a dense, smooth result. Accuracy is also addressed by a low expansion.

- Perfect flowable stone
- Bubble free
- Smooth models

Indication for use: Basing crown and bridge models

SETTING TIME 9-11 min
 WATER/POWDER RATIO 23ml:100g
 COMPRESSIVE STRENGTH 12,500 psi
 EXPANSION 0.08%



BLUE



25LB 1031525
50LB 1031550

BUFF



1031025
1031050

WHITE*



1033025
1033050

Manufactured by Garreco, LLC
430 Hiram Road, Heber Springs, AR 72543, USA
Toll Free: (800) 334-1443
garreco.com

* **WARNING: Cancer** - www.P65Warnings.ca.gov

Garreco, LLC is an ISO 13485 Facility
Garreco is a trademark of Garreco, LLC. Heber Springs, AR. ©2019 All rights reserved.