



Garreco

DIE STONE

VERI-DIE TYPE IV DIE STONE

Produced to handle the many different laboratory techniques that are being used in today's dental labs. With its range in application and versatility, Veri-Die is perfect for almost every requirement. Its physical properties allow Veri-Die to conform to the requirements for each use. Due to its range in particle size, Veri-Die virtually eliminates normal chipping and abrasion. Veri-Die has become the worry free material of choice by providing the technician ample working time and a moderate expansion rate.

- High compressive strength
- Ample working time
- Moderate expansion (0.15%)

Indication for use: Crown and bridge; porcelain; implants; cast partials; high-precision cases

SETTING TIME 10-13 min
 WATER/POWDER RATIO 21ml:100g
 COMPRESSIVE STRENGTH 16,500 psi
 EXPANSION 0.15%



BLUE



GREEN



PINK



IVORY*



120x70g	1401580	1402580	1405080	1403580
25LB	1401525	1402525	1405025	1403525
50LB	1401550	1402550	1405050	1403550

Manufactured by Garreco, LLC
 430 Hiram Road, Heber Springs, AR 72543, USA
 Toll Free: (800) 334-1443
 garreco.com

* WARNING: Cancer - www.P65Warnings.ca.gov

Garreco, LLC is an ISO 13485 Facility
 Garreco is a trademark of Garreco, LLC. Heber Springs, AR. ©2019 All rights reserved.