



VITEX DUPLICATING MATERIAL

Although developed for use in the construction of partial denture frameworks, this versatile hydrocolloid material has proven to be excellent for a variety of dental related tasks. With superior reproduction of detail and endurance for multiple melting cycles. Available in two types (all purpose and regular). All Purpose melts at 190° and holds at 135°. Regular melts at 200° and holds at 125°. Available in Neutral or Blue.

- Versatile
- Superior detail
- Available in 2 colors

Indication for use: Partial denture framework construction

Manufactured by Garreco, LLC
430 Hiram Road, Heber Springs, AR 72543, USA
Toll Free: (800) 334-1443
garreco.com



VITEX DUPLICATING MATERIAL SIZES AVAILABLE

ALL PURPOSE (Melt at 190°F, hold at 135°F)

2 GAL BLUE - 5001015
5 GAL BLUE - 5001050
2 GAL NEUTRAL - 5000015
5 GAL NEUTRAL - 5000050

REGULAR (Melt at 200°F, hold at 125°F)

2 GAL BLUE - 5001515
5 GAL BLUE - 5001550
2 GAL NEUTRAL - 5000115