



## UTILITY WAX



### REGULAR BOXING WAX

A pliable wax that does not need heating yet is easy to work with and is strong and tough. The wax is designed to reduce the possibility of distorting the peripheral roll when boxing an alginate or hydrocolloid impression.

BOXING WAX WHITE - 51510201

BOXING WAX RED - 51510205



### SQUARE ROPES

A soft, pliable, tacky wax that can be used for boxing alginates and hydrocolloids and building up impression trays. It is sufficiently soft to adapt itself to allow the tray to be seated. It will adhere to most any material and can be used for post damming and boxing.

SQUARE ROPES WHITE - 53410000

SQUARE ROPES RED - 53418300



### LARGE UTILITY STRIPS

Large Utility Wax Strips are sticky enough to adhere without softening, but not too sticky to be removed cleanly. These 11" long wax strips are perfect for border molding impression trays and beading impressions to preserve peripheral roll.

LARGE UTILITY STRIPS WHITE - 53710000

LARGE UTILITY STRIPS RED - 53718300



### SMALL UTILITY STRIPS

Small Utility Wax Strips are sticky enough to adhere without softening, but not too sticky to be removed cleanly. These 11" long wax strips are perfect for border molding impression trays and beading impressions to preserve peripheral roll.

SMALL UTILITY STRIPS WHITE - 54710000

SMALL UTILITY STRIPS RED - 54718300



### PERIPHERY STRIPS

Periphery Wax can be used with all types of impression materials for making wax stops, post-damming and lining the periphery of trays.

PERIPHERY STRIPS BLUE - 58711500

PERIPHERY STRIPS WHITE - 58719901