



## RAPIDSIL

### SILICONE DUPLICATING MATERIAL

Lower cost alternative to Accusil with a slightly reduced Shore A hardness of 15 and a faster curing time. Rapidsil offers no expansion or contraction. Due to the use of a platinum base, this material may be poured hundreds of times with equally pleasing results.

- Sold in sets or individually
- 0.0% expansion
- High strength
- Economically-priced

**Indication for use:** Duplicating material

MIXING RATIO..... 1:1 A/B  
WORKING TIME ..... >4 min  
CURING TIME ..... 15-20 min  
ELONGATION OF TEAR ..... ~300%

Manufactured by Garreco, LLC  
430 Hiram Road, Heber Springs, AR 72543, USA  
Toll Free: (800) 334-1443  
garreco.com

Revised on 12/14/20



#### RAPIDSIL SIZES AVAILABLE

1 KG SET (1 ACTIVATOR & 1 BASE) - 6701805  
25 KG SET (1 ACTIVATOR & 1 BASE) - 6701855