



Garreco

SELF CURE

AXCENT ACRYLIC ORTHODONTIC POWDER

Axcent Orthodontic Acrylic has been formulated with added stackability for orthodontic applications. Stock shades may be mixed to make hundreds of unique shades. Micro-fine metallic glitter may be added to any acrylic appliance for added characterization.

- Can be mixed to make unique shades
- Stackable

Indication for use: Orthodontic procedures

CURE TYPE Self Cure
 INJECTABLE No
 MICROWAVABLE No
 POWDER/LIQUID RATIO 2:1



ORTHODONTIC POWDER COLORS AND SIZES AVAILABLE

CLEAR	WHITE	BLACK
4 OZ BOTTLE - 3760104S	4 OZ BOTTLE - 3763004S	4 OZ BOTTLE - 3766304S
4 OZ CASE (4x4 OZ) - 3760104	4 OZ CASE (4x4 OZ) - 3763004	4 OZ CASE (4x4 OZ) - 3766304
1 LB JAR - 3760101S	1 LB JAR - 3763001S	1 LB JAR - 3766301S
1 LB CASE (4x1 LB) - 3760101	1 LB CASE (4x1 LB) - 3763001	1 LB CASE (4x1 LB) - 3766301
5 LB JAR - 3760105S	5 LB JAR - 3763005S	5 LB JAR - 3766305S
5 LB CASE (4x5 LB) - 3760105	5 LB CASE (4x5 LB) - 3763005	5 LB CASE (4x5 LB) - 3766305
20 LB PAIL - 3760120	20 LB PAIL - 3763020	20 LB PAIL - 3766320



MATCHING LIQUID MONOMER

FOR CLEAR ACRYLIC:
Liquid Monomer #4 - Self Cure
 8 OZ BOTTLE - 3860008S
 8 OZ CASE (4x8 OZ) - 3860008 1
 QT BOTTLE - 3860030S
 1 QT CASE - (4x1 QT) - 3860030

FOR COLORED ACRYLIC:
Liquid Monomer #3 - Self Cure
 8 OZ BOTTLE - 3830008S
 8 OZ CASE (4x8 OZ) - 3830008
 1 QT BOTTLE - 3830030S
 1 QT CASE - (4x1 QT) - 3830030

Manufactured by Garreco, LLC
 430 Hiram Road, Heber Springs, AR 72543, USA
 Toll Free: (800) 334-1443
 garreco.com